



3M™ True Definition Scanner

Doctors with a 3M True Definition scanner can send digital impressions directly to Williams Dental Laboratory through the 3M Connection Center.

To send a file:

- 1) Log in to your 3M Connection Center Account that was created when your scanner was installed. If assistance is needed, contact 3M Digital Impression Customer Service at 877-722-528.
- 2) Verify that Williams Dental Laboratory and Williams Dental Laboratory (No Model), are available under the Scanner Laboratory menu.
 - If not, contact the 3M Digital Impression Customer Service (877-722-6528) and request adding the profiles to the scanner menu, or call Williams Dental Laboratory (800-713-5390) ask for Bob.
- 3) Select the Williams Dental Laboratory profile from the drop down menu that matches your case preference:
 - Williams Dental Laboratory (choose this option when requesting models)
 - Williams Dental Laboratory (No Model)* (choose this option when requesting model-less restorations)
- 4) Complete the intraoral scan and select the Williams Dental Laboratory profile that meets your model preference.
- 5) Complete the 3M True Definition online Rx.
- 6) Upload the scan data to the 3M Communication Center.

*Williams Preferred Clients may request a model-less monolithic restoration (zirconium or emax) only by selecting the Williams Dental Laboratory (No-Model) profile that has been previously added to the scanner menu. (see step 2 above).