



Connection Center

For Dental Labs




A digital hub for labs and dentists

Work the way YOU want with digital impressions from the 3M™ True Definition Scanner.

The 3M™ Connection Center is the secure, cloud-based digital hub where dentists upload, store and share prescriptions and digital impressions from the 3M True Definition Scanner.

For traditional work, you simply start with a model made from the digital impression. For CAD/CAM work, you are not boxed into one system. You can rely on trusted and integrated workflows for seamless CAD/CAM production. Or, you can receive an open STL file for import into other systems.

Regardless of how you work, the 3M Connection Center offers the flexibility you need.



3M™ Margin Marking Software is a trusted and integrated digital solution. You can quickly and accurately mark the margin with stereoscopic 3D images.



There's a 3M™ Connection Center package that's right for you.

To assure you get off to a successful start, every dental lab begins as an affiliate lab—which provides you with a working orientation to the 3M™ Connection Center. Once you've completed 30 cases, you can move up to the Standard Lab Package if desired, and have access to the full digital workflow—including the license and training for 3M™ Margin Marking Software. Choose the connection package that's right for your lab.

Affiliate Lab Package

For the dental lab that wants to add a digital option without interrupting the traditional workflow, the Affiliate Package connects you with the basic services you need.

- No up-front investment.
- Confidentially outsource margin marking to a certified margin marking lab (CMML) for a nominal fee. The authorized lab cannot see the name of the dentist or patient.
- Purchase a stereolithography (SLA) model from a trusted model production center. This gives you the ability to create traditional restorations using a digitally produced SLA model.
- Track status of outsourced services with the audit-trail feature of the 3M Connection Center.

Standard Lab Package

The Standard Package provides unlimited access to the 3M Connection Center, and allows you to work with a wide variety of design and production solutions.

- Unlimited access to the 3M Connection Center. This gives you access to scan files, prescription information, and full workflow visibility.
- License for the 3M Margin Marking Software. This gives you the ability to mark margins, set articulation planes, and ditch and die cut models (after completing software training).
- Secure access to open STL files. The 3M Connection Center is the hub for open workflows where you can obtain open STL scan files that can be readily integrated into a variety of digital workflows—allowing you to leverage your existing investment.

Lab Packages	Affiliate	Standard
3M Connection Center lab partner listing	•	•
Trusted SLA model production connection	•	•
Unlimited 3M Connection Center Web access	•	•
Secure file transfer	•	•
3M Margin Marking Software		•
Secure access to open STL files		•
Software upgrades and updates	•	•
Monthly price	FREE	\$79 ^{1,2}

1. One (1) license per 500 cases/annually. Enterprise pricing also available.
2. Billed annually.



Customer Care Center: 1-800-634-2249

www.3M.com/TrueDefinition



3M ESPE

Dental Products

2510 Conway Avenue
St. Paul, MN 55144-1000 USA
1-800-634-2249



70-2013-0558-1
Please recycle. Printed in U.S.A.

3M is a trademark of 3M or 3M Deutschland GmbH. Used under license in Canada. © 3M 2013. All rights reserved.